

# TotalCare® Connect Bed and TotalCare® Bed



## Options and Accessories

**Hill-Rom**

Enhancing Outcomes for Patients and Their Caregivers.™

## Hill-Rom.

### TotalCare® Connect and TotalCare® Bed Options and Accessories

#### Modules

*(Available with the TotalCare SpO<sub>2</sub>RT® surface only)*



#### Percussion & Vibration Module

**product no. P1939E100**

#### Continuous Lateral Rotation Therapy (CLRT) Module

**product no. P1938E100**

#### Low Air Loss Module

*(Available on Rev. L beds and newer only)*

**product no. P1966A**

#### Line Management/IV Poles



#### Line Manager

*(For use on Rev. A-K beds)*

**product no. P1967A**

- Helps organize and manage patient fluid lines




#### Infusion Support System Kit

**product no. K15512**

- The system includes one (P158) ISS pole, two adapters, and one (P159) off-set bar
- The pole will support IV pumps or bags in an upright position
- The head end of the bed has attachment points for two mobile infusion supports

## Line Management/IV Poles *(continued)*



### Infusion Support System

#### product no. P158A


- Consists of a movable and adjustable IV pole that will support IV pumps and bags in an upright position
- Each ISS can support one infusion pump plus two liters of fluid. Safe working load is 20 lbs (9 kg)
- The head end of the bed has attaching points for two mobile ISS



### Set of 2 Permanent IV Poles

#### product no. P1924B110

- Supports up to two IV pumps plus bags
- IV poles attach to the frame near the corners of the headboard
- Can support up to 40lbs (18.1 kg) of total weight
- Retracts and folds down for storage at head of bed



### Removable Telescopic IV Pole (Standard model)

#### product no. P2217

- Installs at any end of the four corners of the bed in the holes provided
- The pole can hold up to 25 lbs (11kg)
- To install, insert pole and rotate in a quarter-turn clockwise; removal is reverse of installation



### Transducer Holder

#### product no. P924

- Use of the Infusion Support System (ISS) is required to use the Transducer Holder
- Transducer safe working load is 2.5 lbs (1.3 kg)



## Hill-Rom.

### TotalCare® Connect and TotalCare® Bed Options and Accessories

## Siderail Pads



### Footboard extenders

*(TotalCare® Connect bed only)*

**product no. P1969A**

### Foot Pad

**product no. P1929A**

- Soft padding with Velcro attachment straps onto foot end of bed
- Provides comfort to patient's feet while resting on footboard

### Siderail Extender

*(For use with Treatment/Short Stay surfaces on the classic TotalCare® bed)*

**product no. P1949EB01**

### Siderail Extender

*(For use with TotalCare SpO<sub>2</sub>RT® surface on classic TotalCare® bed)*

**product no. P1949EB02**



### Siderail Pads with Integrated Extender

*(For use with Treatment/Short Stay surface on the classic TotalCare® bed)*

**product no. P1930EA01**

### Siderail Pads with Integrated Extender

*(For use with TotalCare SpO<sub>2</sub>RT® surface on classic TotalCare® bed)*

**product no. P1930EA02**

- Padded siderail covers (sold only with the Extender) with Velcro to hold in place on the head and intermediate rails
- Cut-away areas for the patient and caregiver control access

## Siderail Pads *(continued)*



### Siderail pads

*(for use with the TotalCare Connect bed - all surfaces)*

- Includes pads for all four siderails

**product no: P1970A**



### Siderail Extender

*(for use with TotalCare Connect bed - any surface)*

**product no:  
P1968A01**

## Traction



### Hill-Rom Patient Helper

#### Patient helper base

**product no. P1842A**

#### Tub Cover for Patient Helper base

**product no. P1848A**


*NOTE: for a complete patient helper assembly, order all 4*

### Headboard for use with Patient Helper

**product no. P1849A**

### Patient Helper Arm with Trapeze *(arm only)*

**product no. P1851A**



### Fracture Frame Adapter Brackets

**product no. P1940A**

- The fracture frame adapter brackets (model 1940) available from Hill-Rom will only work with the fracture frame from Texas Medical Industries (TMI)
- The mount holes for the adapter brackets are located on the weigh frame, under the head section and the thigh section

*NOTE: Fracture Frame Accessories from Zimmer, call: 800-321-5533, ext. 320/324  
Fracture Frame Accessories from TMI, call: 800-332-9556*



Enhancing Outcomes for Patients and Their Caregivers.™

## Hill-Rom.

### TotalCare® Connect and TotalCare® Bed Options and Accessories

## Transport Equipment



### **IntelliDrive® Power Transport** **product no. TC16**

- IntelliDrive® Power Assisted Transport system for TotalCare® bed allows a single caregiver to easily transport the bed in most cases
- Variable speed control
- Available as a factory-installed unit or as a field-upgrade to TotalCare® bed



### **Transport Shelf** *(footboard)* **(factory option)**

- The Transport Shelf provides a platform for small equipment during patient transfer and serves as a writing surface during charting.
- The shelf has a vertical load capacity of 45 lbs (20kg)
- A caregiver can quickly install or remove the shelf in a single step



### **Rack** **product no. P2350A01**

### **Cart** **product no. P2351A01**

### **TotalCare Bed Arm** **product no. P2353A01**

### **Latitude Arm** **product no. P2352A01**

### **Prima Xtend Arm** **product no. P2352A02**

### **Infusion Management Device**

- Facilitates a safe, ergonomic transfer of infusion pumps during patient transfer

## Bed Controls



### Patient Controls

**product no. P1934B (4 button)**

**product no. P1934A01 (11 button)**

- Easy to install using existing siderail hook-up
- Provides head up/down and knee up/down bed functions
- Offers easy patient access to siderail controls



### Nurse Call, Universal TV & Lighting (factory option or field upgrade kit)

- Nurse call, Universal TV and Lighting built into the siderails for easy patient access

## Vertical Oxygen Tank Holder



### Vertical Oxygen Tank Holder

**product no. P27601**

- Attaches to the head end of the articulating frame in a vertical position allowing the tank holders to pivot
- Accommodates one 'E' size oxygen tank with a regulator
- Vertical Oxygen tank holder safe working load is 30 lbs (13.6 kg)



## Hill-Rom.

### TotalCare® Connect and TotalCare® Bed Options and Accessories

#### Accessory Outlet



#### Accessory Receptacle (factory option)

- Receptacle provides up to 12 amps of AC current to provide a convenient source of AC power for non-life support accessory devices
- The receptacle is isolated from the bed system's AC supply and is equipped with an independent power cord per UL code
- Easy access at the foot end of the base frame

#### Casters



#### Caster Options (factory installed)

- 5" Urethane High Mobility Central Locking Caster System (Standard)
- 6" Urethane High Mobility Central Locking Caster System (optional)
- 5" Twin Caster (optional)

#### Seatbelt



#### Seatbelt product no. P1926B

- The seat belt accessory is adjustable
- Holes located on both sides of the bed, accommodate installation of the seat belt

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2009 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

CTG065ra rev 4 3/05/09

**Hill-Rom**

Enhancing Outcomes for Patients and Their Caregivers.™

USA **800-638-2546**  
Canada **800-267-2337**  
[www.hill-rom.com](http://www.hill-rom.com)