

Function, aesthetics and durability
from the partner you trust.



Hill-Rom® Resident® LTC Bed

Hill-Rom

Enhancing Outcomes for Patients and Their Caregivers.™

Enhancing outcomes throughout the care continuum.



People

Clinically-trained sales force



Process

National service network



Technology

Outcome-based solutions

Hill-Rom has been an industry leader and trusted partner in health care since 1929, when our founder, William A. Hillenbrand, decided to “bring the home into the hospital.”

Since then, we have committed ourselves to providing both caregivers and residents high quality, innovative products and solutions, along with support from a dedicated team of clinical experts and service personnel.

Today, our Extended Care Team provides innovative solutions for caregivers and residents through our clinically-trained sales force and a national network of Hill-Rom service technicians. We offer a full portfolio of post-acute products, including bed frames, prevention and therapy surfaces, bariatric solutions, furniture and clinical support and education services.

You can count on the people, processes and technology of Hill-Rom—in hospitals, extended care facilities and in the home.

“Hill-Rom has been an industry leader and trusted partner in health care since, when our founder, William A. Hillenbrand, decided to ‘bring the home into the hospital.’”

Hill-Rom offers a wide array of products and services to complement the Hill-Rom® Resident® LTC Bed, including:

Prevention mattresses	Therapy surfaces	Furniture
<ul style="list-style-type: none"> • NP100 • NP50 	<ul style="list-style-type: none"> • Synergy® Air Elite • P300 	<ul style="list-style-type: none"> • Overbed tables • Bedside cabinets



Technical Specifications – Hill-Rom® Resident® LTC Bed

Feature	Dimension
Width – Siderails Up	42" (106.7 cm)
Width – Siderails Down	40.25" (91.4 cm)
Length – Frame Only	86.5" (219.7 cm)
Length – With Bumpers in Low Position	91.5" (232.4 cm)
Feature	Dimension
Deck Height Range	13.875 - 30" (35.2 - 76.2 cm)
Approximate Bed Weight	275 lbs (120.2 kg)
Sleep Deck	36" x 80" (91.4 x 203.2 cm)

Hill-Rom® Resident® LTC Bed

The Resident® LTC Bed was designed specifically for long term care. It provides a pleasing combination of function, aesthetics, durability, and value to enhance the well-being of both the resident and caregiver. Several options are available to meet your specific needs. Please contact your local Hill-Rom representative for additional information regarding these options.



Resident LTC Bed with optional Turtlecreek™ head and footboards and optional assist bars



The standard head and knee elevation with one-touch Auto Contour enhances caregiver safety and allows for easier patient positioning.

Additional standard features:



New

Head of Bed angle indicators help ensure proper elevation.



Footend control lockout panel helps create a safer resident environment.



Locking casters allow for easy relocation and positioning within the room.



Large, easy-to-see, easy-to-operate siderail-mounted controls enhance resident comfort and independence.



Electric high/low movement facilitates resident egress and provides comfortable working height for caregivers.



Standard drainage hooks and IV pole mounts enhance safety for both caregivers and residents.

The Hill-Rom® Resident® LTC Bed accommodates a wide variety of prevention mattresses and treatment surfaces.

Head/Footboard Styles



New

Coventry™



New

Turtlecreek™



New

Jamesburg™

The Coventry™, Turtlecreek™ and Jamesburg™ styles offer attractive styling, durable, thermoformed construction, and attractive colors that match the Hill-Rom® Art of Care® furniture line.

Head/Footboard Color Options (Available on Coventry™, Turtlecreek™ and Jamesburg™ Styles)



American Natural

Red Oak

Solara Oak

Heywood Cherry

Camby Cherry

Newport Apple

Acajou Mahogany

Clove



Century



Hearthside

Head/Footboard Color Options (Available on Century and Hearthside Styles)



Natural Maple

Honey Maple

Light Oak

Medium Oak

Wild Cherry

Shaker Cherry

Dark Cherry

Montana Walnut

Package Configurations

	LTC320	LTC310	LTC300	LTC250	LTC205	LTC105	LTC100
Head rails with controls	•	•	•				
Head rails with no controls				•			•
Foot rails with no controls	•						
Assist Bar (foot end)		•					
Two floor breaks						•	•
Four floor breaks	•	•	•	•	•		
Pendant				•	•	•	•
No head/foot rails					•	•	

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2009 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.



Enhancing Outcomes for Patients and Their Caregivers.™

USA 800-445-3730
Canada 800-267-2337

www.hill-rom.com