

# Hill-Rom®

## 1000 Bed



## Introducing the Hill-Rom® 1000 Bed.

Value minded doesn't mean sacrificing innovation and style. The Hill-Rom 1000 Bed fits the need for both patient and caregiver — dedicated to safety, comfort and efficiency.

### At a glance

- One-button Dining Chair® position for patient comfort and safety
- Patient pendant is stowable in any siderail
- Shearless Pivot®, mechanism with auto contour means less need for caregivers to reposition patients
- Higher headboard makes moving the bed easier
- Easy-to-read Line-of-Site® angle indicators for degree of head section and Trend and Reverse Trend
- *Instant CPR* release lever includes a dampened descent, gently lowering the patient and reducing caregiver back strain
- Point-of-Care® controls save caregiver steps
- Four dual-locking casters keep the bed in place
- Bed exit system
- Bed and mattress are designed to address main entrapment areas
- Optional built-in scale eliminates the need to move patients to a separate scale

**Hill-Rom®**

Enhancing Outcomes for Patients and Their Caregivers.™

**Hill-Rom®**

1000 Bed

## Standard Features Unique to Hill-Rom

### Patient Safety

Your patients' safety and comfort are your primary concern. The Hill-Rom 1000 Bed incorporates several unique features to help you protect and care for your patients.



#### **One-button Dining Chair. feature**

- One-button control moves patient in and out of chair position.
- The bed maintains a low, safe height while putting patient in a chair position.



#### **Stowable patient pendant**

For ease of use, patients have full access to bed controls regardless of rail position.



#### **Unique Brake Bar at the Footend of the Bed**

- The brake bar is easily accessible even when siderails are down.
- The brake bar is located away from patient egress points.



#### **Shearless Pivot. feature**

This unique feature can reduce the need to reposition patients with the touch of a button.

### Caregiver Ergonomics

With features that help provide safer patient mobility and ease of use, the Hill-Rom 1000 Bed was designed with caregiver safety as a priority.



#### **Non-moving headboard**

- Headboard remains at ergonomic height during all bed articulations.
- Stationary headboard height means no more concerns about IV poles breaking the lights.

#### **Line of Sight. angle indicator**

The indicator is located on the siderail for easy visibility.

**Hill-Rom®**

1000 Bed

## Caregiver Convenience

The Hill-Rom 1000 Bed controls are created to be easy-to-use and promote efficient patient care.



### **Point-of-Care® controls**

Identical, easy-to-read caregiver controls, including scale, are located on both sides of the bed.

### **Battery back-up**

- There is no need to find an outlet when repositioning the bed during patient transfer.
- Battery back-up automatically charges when the bed is plugged in.

## Additional Hill-Rom 1000 Bed Advantages

The Hill-Rom 1000 Bed meets FDA Guidance on Hospital Bed System Dimensions (HBSW) to reduce the potential for siderail entrapment.

### **Scale**

- Built-in scale eliminates the need for sling scale procedures helping to improve patient safety.
- Located on the siderail for convenient, easy-to-read controls.
- Scale can weigh in any position.

### **Bed Exit System**

- Setting the bed exit system, combined with falls prevention protocols, can help to reduce patient falls.
- Newer, more reliable technology means fewer false calls.
- This system is single-mode with both tone and volume controls.



## Technical Specifications

### Hill-Rom 1000 Bed

<b>Feature</b>	<b>Dimension</b>
Frame width.....	39.74" (100.95 cm)
Frame length with head and footboard (w/ bumpers) normal.....	90.04" (228.71 cm)
Sleep deck width.....	36.00" (91.44 cm)
normal length.....	80.00" (203.20 cm)
Bed clearance under frame.....	5.00" (12.70 cm)
Optional Scale System accuracy.....	+/- 2.2 lbs or 1.1% of patient weight w/e is greater
weigh in any position.....	yes
scale repeatability.....	2% of patient weight

<b>Feature</b>	<b>Dimension</b>
Low position floor to deck.....	15.75" (42.55 cm)
High position floor to deck.....	32.5" (82.55 cm)
Max head elevation.....	65 degrees
Max Trend and reverse Trend.....	16 degrees
Caster size.....	6" (15 cm)

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2008 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

153249 rev 1 4/21/08



Enhancing Outcomes for Patients and Their Caregivers.™

USA 800-445-3730  
Canada 800-267-2337

[www.hill-rom.com](http://www.hill-rom.com)